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
The basic principle

One of the most important aspects in the isolation of RNA is to prevent degradation during the isolation procedure. With the **NucleoSpin® RNA** methods, cells are lysed by incubation in a solution containing large amounts of chaotropic ions. The RA1 lysis buffer immediately inactivates RNases – which are present in virtually all biological materials – and creates appropriate binding conditions which favor adsorption of RNA to the silica membrane. Contaminating DNA, which is also bound to the silica membrane, is removed by an rDNase solution which is directly applied onto the silica membrane during the preparation (RNase-free rDNase is supplied with the kit). Simple washing steps with two different buffers remove salts, metabolites and macromolecular cellular components. Pure RNA is finally eluted under low ionic strength conditions with RNase free H₂O (supplied).

Component	Hazard contents	GHS symbol	Hazard phrases	Precaution phrases
Inhalt	Gefahrstoff	GHS-Symbol	H-Sätze	P-Sätze
rDNase, RNase-free	rDNase, lyophilized <i>rDNase, lyophilisiert</i> CAS 9003-98-9	 DANGER GEFAHR	317, 334	261, 280, 302+352, 304+340, 333+313, 342+311, 363
RA1	guanidinium thiocyanate 30–60 % <i>Guanidinthiocyanat 30–60 %</i> CAS 593-84-0	 WARNING ACHTUNG	302, 412, EUH031	260, 273, 301+312, 330

1	Homogenize sample		30 mg		
2	Lyse cells		350 µL RA1 3.5 µL β-mercaptoethanol Mix		
3	Filtrate lysate		11,000 x g, 1 min		
4	Adjust RNA binding conditions		350 µL 70 % ethanol Mix		
5	Bind RNA		Load sample 11,000 x g, 30 s		
6	Desalt silica membrane		350 µL MDB 11,000 x g, 1 min		
7	Digest DNA		95 µL DNase reaction mixture RT, 15 min		
8	Wash and dry silica membrane		1 st wash 200 µL RAW2 2 nd wash 600 µL RA3 3 rd wash 250 µL RA3 1 st and 2 nd  11,000 x g, 30 s 3 rd  11,000 x g, 2 min		
9	Elute highly pure RNA		60 µL RNase-free H ₂ O 11,000 x g, 1 min		

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pGEM®-T Vector Systems

The pGEM®-T Vector Systems are convenient systems to clone PCR products. The pGEM®-T Vector is prepared by cutting the pGEM®-5Zf(+) Vector with EcoRV and adding a 3' terminal thymidine to both ends. These single 3'-T overhangs at the insertion site greatly improve the efficiency of ligation of a PCR product into the plasmid by preventing recircularization of the vector and providing a compatible overhang for ligation of PCR products generated by thermostable polymerases that add a single deoxyadenosine, in a template-independent fashion, to the 3'-ends of amplified fragments.

Features - Benefits

- **Rapid Ligation:** The provided 2X Rapid Ligation Buffer allows reactions to be completed in 1 hour at room temperature.
- **Blue/White Screening:** T7 and SP6 RNA polymerase promoters flank a multiple cloning region within the α -peptide coding region for β -galactosidase. Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates.



Table 1. Genotypes of *E. coli* Competent Cells Offered by Promega.

Strain	Genotype
BL21(DE3)pLysS	F ⁻ , <i>ompT</i> , <i>hdsS_B</i> (<i>r_B</i> ⁻ , <i>m_B</i> ⁻), <i>dcm</i> , <i>gal</i> , λ (DE3), pLysS, Cm ^r
HB101	F ⁻ , <i>thi-1</i> , <i>hdsS20</i> (<i>r_B</i> ⁻ , <i>m_B</i> ⁻), <i>supE44</i> , <i>recA13</i> , <i>ara-14</i> , <i>leuB6</i> , <i>proA2</i> , <i>lacY1</i> , <i>galK2</i> , <i>rpsL20</i> (<i>str^r</i>), <i>xyl-5</i> , <i>mtl-1</i>
JM109	<i>endA1</i> , <i>recA1</i> , <i>gyrA96</i> , <i>thi</i> , <i>hdsR17</i> (<i>r_k</i> ⁻ , <i>m_k</i> ⁺), <i>relA1</i> , <i>supE44</i> , Δ (<i>lac-proAB</i>), [F ['] , <i>traD36</i> , <i>proAB</i> , <i>lacI^q</i> , <i>lacZΔM15</i>]
Pro 5-alpha	<i>fhuA2</i> , Δ (<i>argF-lacZ</i>), U169, <i>phoA</i> , <i>glnV44</i> , ϕ 80, Δ (<i>lacZ</i>)M15, <i>gyrA96</i> , <i>recA1</i> , <i>relA1</i> , <i>endA1</i> , <i>thi-1</i> , <i>hdsR17</i>

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
5A : Molecular model of T₃ target genes activation

Gene regulation by T₃ hormone is achieved through activation of the Thyroid Hormone receptor (TR α 1). Heterodimeric RXR/TR α 1 complex can bind to thyroid hormone response elements (TRE) in DNA. In the presence of T₃ hormone, a conformational change in TR α 1 occurs permitting activation of target genes.

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